

SLC25A27-solute carrier family 25, member 27 Gene

Reference: PTXBC033091

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UCP4; mitochondrial uncoupling protein 4; UCP 4; uncoupling protein 4; solute carrier family 25 member 27		

For research use only.

GENE DESCRIPTION	solute carrier family 25, member 27
GENE ID	9481
ORIGIN SPECIES	Human
ACCESSIONS	BC033091
VECTOR	pDONR223
NCBI SYMBOL	SLC25A27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgtccggaggaggaggaggtctttgcccgtgaccagagatggccccgagcgagcaaatctactgtccggtgcgcggctacgtggcggagctagcaaccttccccctggatctcacaactcgactccaaatgcaaggagaagcagctctgtctcgggtgggaga
cgggtgcaagagaatctgccccctatagggaatgggtgcgcacagccctaggatcatgaagaggaaggctttctaaagcttggcaaggagtgacacccgccatttacagacacgtagtgtattctggaggtcgaatgggtcacatatgaacatctccgagaggtgtgttggcaaaagt
aagatgagcattatcccccttggaaatcagtcattggaggatgatggctggtgtattggccagtttttagtcaatccaactgacctagtgaagggttcagatgcaaatggaagaaaaagaaactggaaggaaaaccattgcgatttcgtgtgtacatcatgcatcttgcataatcttagctg
aaggaggataacgagggcttggcaggctgggtaccataatacaaaagacagcactggtgaatatgggagatttaaccactatgatacagtgaaacactactgtgtattgaatacaccactgaggacaatatcatgactacgggtttcaagtgtctgtcggatctcacaaggccat
ccaatga