

## HBEGF-heparin-binding EGF-like growth factor Gene

Reference: PTXBC033097

|                |                                                                                                                                                                                                                                                                |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | DTR; DTS; DTSF; HEGFL; proheparin-binding EGF-like growth factor; diphtheria toxin receptor (heparin-binding EGF-like growth factor); diphtheria toxin receptor (heparin-binding epidermal growth factor-like growth factor); heparin-binding epidermal growth |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | heparin-binding EGF-like growth factor                                                            |
| <b>GENE ID</b>          | 1839                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC033097                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | HBEGF                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaagctgctgccgtcggtggtgctgaagctcttctgctgcaagctctctcgccactggtgactggcgagagcctggagcggcttcggagagggctagctgctggaaccagcaacccggaccctccactgtatccacggaccagctgctaccctaggaggcgccgggaccgga  
aagtcctgactgcaagaggcagatctggaccttttgagagtcactttatctccaagccacaagcactggccacaccaacaaggagagcacgggaaaaaagaagaagaaggcaaggggctagggaagaaggagggaccatgtcttcggaaatacaggacttctgcatcctg  
gagaatgcaaatatgtgaaggagctccgggctccctcctgcatctgccaccgggttaccatggagagaggtgtcatgggtgagcctccagtgaaaaatcgcttatatactatgaccacacaacacatcctggccgtggtgctgtggtgctgtcatctgtctgtgctgctggtgctggtg  
ggcttctcatgtttaggtaccataggagagaggttatgtgtgaaaaatgaagagaagtgaagtggcgatgactaattccactga