

hCG_17241-flamingo-like Gene

Reference: PTXBC033112

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	flamingo-like
GENE ID	728032
ORIGIN SPECIES	Human
ACCESSIONS	BC033112
VECTOR	pDONR223
NCBI SYMBOL	hCG_17241
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgagctcttagggcagcaaaagccttcagggtgatgtccctcgggtccacagccaggcctccgaggctgagcactgagcgtgtgtgtgtgaccagaacatcaccactgggcagtgccctcctggccccgccaccagggcagcgtgggagcagatccagcggagc
gagggtggcagggcagcagctgtccggcgtcagggtctactcagcaatgtggcagcaacgtgcagtggaactacctgcagccctttgtcatgtcaccaccaacatgagtaaggcactggcgtgctgtgggggtggagcctgcaaccctcaggggtgcgccaggagctgggt
ccccatgggtggagcaggtgggtcccgcccccaactggccatgctcttagaaaaatga