

TMEM199-transmembrane protein 199 Gene

Reference: PTXBC033113

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf32; CDG2P; VMA12; VPH2; transmembrane protein 199		
For research use only.			

GENE DESCRIPTION	transmembrane protein 199
GENE ID	147007
ORIGIN SPECIES	Human
ACCESSIONS	BC033113
VECTOR	pDONR223
NCBI SYMBOL	TMEM199
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctcctttgttcggcgagcgattggtgctgtgctttgggccccggcgaggagcagagcggtaccgccgaaagctgcggcgccgagcttgaggccgctggggaagaagcacaaggcggtgatagctccagtgcccccaacgcttggtttcttcgtctca
tccgggatctgcaccagcatctgagagaaaggattccaaactatccatgagctcctagaaggcagtgaaatctatctccagaggtgtgaagcctccacggaaccagaactagtgcggcgctggagaagattaagatacagctggccaatgaggaaataaacggatcaccc
gcaacgtcactgtcaggatacaagacatggtgggactctcagcgacctgggaaagcaagtgaagcattgaaggctctggtcatcaccatctcaatttcattgtcacggtggtgtgctgccttcgtctgcacttaccttggagccaatatacttcacagaaatggcctcgccggtgctagctg
cattgatcgtcgctctgtggtgggtctgcccagctgtatgtcatggtgcgggcaatggaaggcgagctgggagaactgtaa