



Reference: PTXBC033199

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 616		
For research use only.			

GENE DESCRIPTION	zinc finger protein 616
GENE ID	90317
ORIGIN SPECIES	Human
ACCESSIONS	BC033199
VECTOR	pDONR223
NCBI SYMBOL	ZNF616
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr