



Reference: PTXBC033211

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MG53; tripartite motif-containing protein 72; mitsugumin-53; tripartite motif containing 72, E3 ubiquitin protein ligase; tripartite motif containing 72		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 72
GENE ID	493829
ORIGIN SPECIES	Human
ACCESSIONS	BC033211
VECTOR	pDONR223
NCBI SYMBOL	TRIM72
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgtcggcctgcgccggcctcctgcaccaggagctgctgcgcgcctgctgcctgcagctgttcgacgcgccgcgcagaccgagtcgcggccacagtttctgcgcgcctcctgagccgcgtgcggccgggagccggcgcggatggccaccgttctctcccctctgctccaggccccc
acgcggcccgacggcactcagcaccaactcgcagctgcgcgcgcctgttgagggggcgtgccagggtgccgagggccactgcgaggagcacctgcgaccgcctgagcatctactgcgagcaggaccgcgcgtgctgtgcggagtgtgcgcctcactcggctcgcaccgcggtcat
cgcctctcctcctccgcggagcccgacgcacgcctcaagacacagctgccacagcagaactcagctgcagtcgagcgcagtcgctgctgaaggagaagatgtggcgtctgagcatcagctgttgagagtggagagacagtgcctcagttccggggggccgtggggagcagc
ctgggcaagatcctgctgtctcgtcgtcactgggggcctctggaccgcgcagagcgtgtacggctgagcgcagctgccttgcgtccgcggagctgcggcgtgaactcttacctggagcagctgcgcagatggagaagctctgcggagcgtgtgcggacaagccg

Page 1/1