



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SIRPD-signal-regulatory protein delta Gene

Reference: PTXBC033502

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

PTPNS1L2; dJ576H24.4; signal-regulatory protein delta; SIRP-delta; protein tyrosine phosphatase, non-receptor type substrate 1-like 2; signal regulatory protein delta

For research use only.

GENE DESCRIPTION	signal-regulatory protein delta
GENE ID	128646
ORIGIN SPECIES	Human
ACCESSIONS	BC033502
VECTOR	pDONR223
NCBI SYMBOL	SIRPD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]