

CXorf41-chromosome X open reading frame 41 Gene

Reference: PTXBC033510

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf41; NYSAR97; protein PIH1D3; PIH1 domain-containing protein 3; sarcoma antigen NY-SAR-97; PIH1 domain containing 3		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 41
GENE ID	139212
ORIGIN SPECIES	Human
ACCESSIONS	BC033510
VECTOR	pDONR223
NCBI SYMBOL	CXorf41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaatctgaaaatatggattctgaaaatatgaagacagaaaatatggaatcctcaaatgtagacttgagagtggtttctcagttacagctctggaagccctctaaagctacttaactctgaagaagaggatgattctgactatggacagacaaatggttatctactattggagccatgggtcct
gggaatattggaccacccaatafagaagagctcaaaagtcacccctgaaaccagcagagaaaataatgaggacatctggaattcagaagagattccagaagagcagaatatgatgatattgggatgttagagaaatcccagagtagagattatattcagacagcaggtgggaactgaa
gatataatttttaggggtgtcaaaaaggactcctcaacaggtgtgtgagtgactagtggctaaaaataaattgccaaatacaaacctctgatattcaaatgatccaggaacaatcctgacctctgactcctcagaagaagctgttgataactctcctgagctgtggaatgtaccagt
ccaaagcattctatataccagagactgaaactcttgaatcactatgactatgaaaagagagtttagatttcttctga