



Reference: PTXBC033571

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR92; GPR93; KPG_010; LPA5; lysophosphatidic acid receptor 5; G protein-coupled receptor 92; G-protein coupled receptor 93; LPA receptor 5; LPA-5		

For research use only.

GENE DESCRIPTION	lysophosphatidic acid receptor 5
GENE ID	57121
ORIGIN SPECIES	Human
ACCESSIONS	BC033571
VECTOR	pDONR223
NCBI SYMBOL	LPAR5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1