



Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

Reference: PTXBC033647

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	3M1; CUL-7; KIAA0076; dJ20C7.5; cullin-7; cullin 7		
For research use only.			

GENE DESCRIPTION	cullin 7
GENE ID	9820
ORIGIN SPECIES	Human
ACCESSIONS	BC033647
VECTOR	pDONR223
NCBI SYMBOL	CUL7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr