

BOLL-bol, boule-like (Drosophila) Gene

Reference: PTXBC033674

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein boule-like; bol, boule-like; boule-like RNA binding protein; boule homolog, RNA binding protein		

For research use only.

GENE DESCRIPTION	bol, boule-like (Drosophila)
GENE ID	66037
ORIGIN SPECIES	Human
ACCESSIONS	BC033674
VECTOR	pDONR223
NCBI SYMBOL	BOLL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaaacagattcattatctccatccctaatctgtgtcacctgtgctttgaataaccaacaagtgccccagatatggaacagtgatccctaatcgatctttgtaggagggaattgatttaagacaaacgaaagtgatttaagaaaatttttccagtagtggtctgtgaaagaagtgaag
attgtaaatgacagagctggagatccaaagggtatggttcgtcacttttgaacacaaagaagatgcacaaaaattttacaagagctgaaaaacttaattataaggataagaagctgaacattggctcagcaataagaaaacaacaagtagggatccctgttctagtataatgccagcag
ctggacaatgtatctaacaactcaactggatatttatactaccataatggtgtgttttctactccagaggttaacttcggccaccgcttggccttcacgttctgtatgtagctccccgtgatggttagctcagccatttatcagcaacctgcatactaccaggccaccacaca
gtatttaccaggacagtgccagtggtgttcctcagccttctgctcctctctgtcctcattctatactgcacacctctgaggttattatcaaccagtggaaattgcacaggatggtggatgtgttctctccactgtctctgatgaaacttcagttccagagccttattctgatcatggagttcaag
caacatatcaccaggtttatgctccaagtgccatcactatgcctgcgctgtgatgcagcctgagccaattaaaacagtggtggagcattcattataa