



Reference: PTXBC033691

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYP11S1; cytochrome P450 2S1; cytochrome P450, family 2, subfamily S, polypeptide 1; cytochrome P540, subfamily IIS, polypeptide 1; cytochrome P450 family 2 subfamily S member 1		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 2, subfamily S, polypeptide 1
GENE ID	29785
ORIGIN SPECIES	Human
ACCESSIONS	BC033691
VECTOR	pDONR223
NCBI SYMBOL	CYP2S1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]