

## BCL2L11-BCL2-like 11 (apoptosis facilitator) Gene

Reference: PTXBC033694

|                |                                                                                                                                                                                                   |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | BAM; BOD; bcl-2-like protein 11; BCL2-like 11 (apoptosis facilitator); bcl-2 interacting mediator of cell death; bcl-2 interacting protein Bim; bcl-2-related ovarian death agonist; BCL2 like 11 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | BCL2-like 11 (apoptosis facilitator)                                                              |
| <b>GENE ID</b>          | 10018                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC033694                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | BCL2L11                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcaagcaacctctgatgtaagttctgagtgaccgagaaggttagacaattgcagcctgcggagaggcctccccagctcagacctggggccctacctccacagacagagccacaaggtaacctgaaggcaatcacggaggtgaaggggacagctgccccacggcagc  
cctcaggggcccgctggccacacctgccagccctggccctttgctaccagatccccgctttcatctttagagaagatcctccctgctgtctgatccctcagtggtatttctctttgacacagacaggagcccgacccatgagttgtgacaaatcaacacaaacccaagtctctcttc  
caggcctcaaccactatctcagtgcaatggctccatgaggcaggtgaacctgcagatatgcgccagagatatggatcgcccaagagttgcggcgtattggagacgagtttaacgcttactatgcaaggagggtattttgaataattaccaagcagccgaagaccaccacgaatggtt  
atcttacgactgttacgttacattgtccgctggtgtggagaatgcattga