



Reference: PTXBC033740

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDS3; HXBL; TENX; TN-X; TNX; TNXB1; TNXB2; TNXBS; VUR8; XBS; tenascin-X; growth-inhibiting protein 45; hexabrachion-like protein; tenascin XB1; tenascin XB2; tenascin XB		

For research use only.

GENE DESCRIPTION	tenascin XB
GENE ID	7148
ORIGIN SPECIES	Human
ACCESSIONS	BC033740
VECTOR	pDONR223
NCBI SYMBOL	TNXB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]