

NDUFV3-NADH dehydrogenase (ubiquinone) flavoprotein 3, 10kDa Gene

Reference: PTXBC033766

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-10k; CI-9kD; NADH dehydrogenase [ubiquinone] flavoprotein 3, mitochondrial; NADH dehydrogenase (ubiquinone) flavoprotein 3, 10kDa; NADH-ubiquinone oxidoreductase 9 kD subunit; NADH-ubiquinone oxidoreductase 9 kDa subunit; NADH-ubiquinone oxidoreductase		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) flavoprotein 3, 10kDa
GENE ID	4731
ORIGIN SPECIES	Human
ACCESSIONS	BC033766
VECTOR	pDONR223
NCBI SYMBOL	NDUFV3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgccccgtgtttgctgcggcaaggacgagccggggcgtgaagactatgctccaggagccaggtgtttcaggacttgcctctacggtttcttctgctgcggcaatcagggaagagtgaaaaggggtcagccacagaattccaagaagcaagtcacacaaaaagccagccccc
agtgccctgctgagccgtttgacaacactacctacaagaacctgcagcatcatgactacagcagctacacctcttagacctcaacctgaacctcaaaatcaggatgcctcagccctctcagccgggagtcacctcgacactga