

TPD52L3-tumor protein D52-like 3 Gene

Reference: PTXBC033792

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NYDSP25; tumor protein D55; protein kinase NYD-SP25; testis development protein NYD-SP25; testis tissue sperm-binding protein Li 87P; tumor protein D52 like 3		

For research use only.

GENE DESCRIPTION	tumor protein D52-like 3
GENE ID	89882
ORIGIN SPECIES	Human
ACCESSIONS	BC033792
VECTOR	pDONR223
NCBI SYMBOL	TPD52L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccacatgccaggacagacacctctgtgggcacatatgaatccactcgaacttgaggatctgacagagcccagcaagagagctcaaaacaaactcactaaattggaggctgaaattgtaacctacgccactactggcagccaaagagagacgtgtggggaact
caagaggaagttaggcctcaccgccttggtaggctgagacagaatctgccaagagctgcttgatgttcagggtctccaacacctatgtgaacagaagacatcagctgctgtgccatgggcactctcatctgcaggaagcttgaggcgtgaagaagtcggccacactcagatctt
ttgaaggattgatctcaataatacacgttaaatcaaggaggaattaa