

HDHC3-HD domain containing 3 Gene

Reference: PTXBC033794

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	(ppGpp)ase; MESH1; guanosine-3',5'-bis(diphosphate) 3'-pyrophosphohydrolase MESH1; HD domain-containing protein 3; metazoan SpoT homolog 1; penta-phosphate guanosine-3'-pyrophosphohydrolase; HD domain containing 3		

For research use only.

GENE DESCRIPTION	HD domain containing 3
GENE ID	374659
ORIGIN SPECIES	Human
ACCESSIONS	BC033794
VECTOR	pDONR223
NCBI SYMBOL	HDHC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggctctgagggcgccagctgctggaggtgccgacttcgggctcgaagcaccggcagcagcggcggaaggaccccgaggggacccctacatcaaccaacccatcgggtggcacggatcctgaccacagagcggggaatcactgacattgtggtgttacaggcggccc
tgcctcatgacacgggtggagacacagacacacccctggatgaggtggagctacatttggggcacaagtgcggcgccctgggtggagggtaacagatgacaagacttgcgaagctggagagaaagagctgcaggtggagcaagcggccacagtagccccggggccaaac
tggatgaagctggcagacaagctgtacaatctgagggacctgaatcgctgcacccagaggtaaaaatacaatag