



Reference: PTXBC033801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIG10; EJM7; GEFSF5; gamma-aminobutyric acid receptor subunit delta; GABA(A) receptor subunit delta; GABA(A) receptor, delta; GABA-A receptor, delta polypeptide; gamma-aminobutyric acid (GABA) A receptor, delta; gamma-aminobutyric acid type A receptor del		

For research use only.

GENE DESCRIPTION	gamma-aminobutyric acid (GABA) A receptor, delta
GENE ID	2563
ORIGIN SPECIES	Human
ACCESSIONS	BC033801
VECTOR	pDONR223
NCBI SYMBOL	GABRD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1