

ABCC9-ATP-binding cassette, sub-family C (CFTR/MRP), member 9 Gene

Reference: PTXBC033804

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABC37; ATFB12; CANTU; CMD1O; SUR2; ATP-binding cassette sub-family C member 9; ATP-binding cassette transporter sub-family C member 9; ATP-binding cassette, sub-family C (CFTR/MRP), member 9; sulfonylurea receptor 2; ATP binding cassette subfamily C membe		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family C (CFTR/MRP), member 9
GENE ID	10060
ORIGIN SPECIES	Human
ACCESSIONS	BC033804
VECTOR	pDONR223
NCBI SYMBOL	ABCC9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccttcatttggtaacaacattctcatataatcaacgatgggtgactacaaaattcctgcttggatgccctcaacctggccctcatgtcttctgtttatcaactttccaatattgtttatgggtggggagccaaagctcaaaagtacaaattcaccacaacacatggettcatttccgggacataacctgagatggattcttattcgtctcctgtttgtgcatgtctgtgaaatagcagaaggcattgtttcagactcgcggcggaatcaaggcacctccacctctttatgccagccgtgatgggattcgttgcactacaacatcgatagtgattatcataatcgaacatcaaatctctaaattacttttaggttttagccatcatgagagtaaaagctgctatcaacattcttaa