



Reference: PTXBC033805

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NLT; OAT2; solute carrier family 22 member 7; hOAT2; liver-specific transporter; novel liver transporter; organic anion transporter 2; solute carrier family 22 (organic anion transporter), member 7		

For research use only.

GENE DESCRIPTION	solute carrier family 22 (organic anion transporter), member 7
GENE ID	10864
ORIGIN SPECIES	Human
ACCESSIONS	BC033805
VECTOR	pDONR223
NCBI SYMBOL	SLC22A7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr