

TRIM73-tripartite motif-containing 73 Gene

Reference: PTXBC033812

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TRIM50B; tripartite motif-containing protein 73; tripartite motif-containing protein 50B; tripartite motif containing 73		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 73
GENE ID	375593
ORIGIN SPECIES	Human
ACCESSIONS	BC033812
VECTOR	pDONR223
NCBI SYMBOL	TRIM73
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccgagaccctaccgacggcgccatctctgtgggcccacgagggcactcatcgtttcgttccaaccgcccgaacggccctggccctcgcgagctctggaactacagaggtcgacgggtgagttgccaggtgtggcccgtaacggagcgcacaaaacatgatggga
cacgtaacgggaccacacagggcacattgggcacttgaggggcgcgaggtggcggcacgtaa