



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CADM3-cell adhesion molecule 3 Gene

Reference: PTXBC033819

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

BigR; IGSF4B; NECL1; Necl-1; synCAM3; cell adhesion molecule 3; TSLC1-like 1; brain immunoglobulin receptor; dendritic cell nectin-like protein 1 short isoform; immunoglobulin superfamily member 4B; synaptic cell adhesion molecule 3

For research use only.

GENE DESCRIPTION	cell adhesion molecule 3
GENE ID	57863
ORIGIN SPECIES	Human
ACCESSIONS	BC033819
VECTOR	pDONR223
NCBI SYMBOL	CADM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr