

DNAJC3-DnaJ (Hsp40) homolog, subfamily C, member 3 Gene

Reference: PTXBC033823

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACPHD; ERdj6; HP58; P58; P58IPK; PRKRI; dnaJ homolog subfamily C member 3; DnaJ (Hsp40) homolog, subfamily C, member 3; ER-resident protein ERdj6; endoplasmic reticulum DNA J domain-containing protein 6; interferon-induced, double-stranded RNA-activated p		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 3
GENE ID	5611
ORIGIN SPECIES	Human
ACCESSIONS	BC033823
VECTOR	pDONR223
NCBI SYMBOL	DNAJC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggccccggctccgtgaccagccggctgggctcgttccctctcgtagtcctggtgatctgcagtagcgaaggctgctgaatgtggagtaaatgcagatgttgagaacatctgaattgggcaagaataacttcagctggacagctagctgatcttctcagtttcagct
gccgtagatgggtgacctgataactatattgcttattatcgagggtcactgtcttttagctatggcgaatcaaaagctgcacttctgatttaactaaagtgaattcaattgaagatggacttcactgcagcaagattacagagaggtcacttattactcaacaaggaaaacttgatgaagcaga
agatgattttaaaaaatgggtgttctctctctgttgggactccagcgttctctcctagatgatctatactgtctattttggttcttctatgaagaaggtgacttcagatgcctctcgtcagccatctctgaatgccttctcagctcttaatttgatgaagttcaattgttgatttcttttacttc
tgtggactgtgcgtttggtgtcatgtctaagaagtattcactacgccgtaggatctaaaactttctcatatctctaaaagtttatggtttgtttttatttttaagccatagctatttgagttaa