

CHMP4B-chromatin modifying protein 4B Gene

Reference: PTXBC033859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf178; CHMP4A; CTPP3; CTRCT31; SNF7; SNF7-2; Shax1; VPS32B; Vps32-2; dJ553F4.4; charged multivesicular body protein 4b; SNF7 homolog associated with Alix 1; Snf7 homolog associated with Alix 1; chromatin modifying protein 4B; chromatin-modifying pro		

For research use only.

GENE DESCRIPTION	chromatin modifying protein 4B
GENE ID	128866
ORIGIN SPECIES	Human
ACCESSIONS	BC033859
VECTOR	pDONR223
NCBI SYMBOL	CHMP4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtgttcgggaagctgttcgggctgaggggtaagccggcaaggggcgccgaccccccaggagggccatccagcggtcggggacacggaagagatgtaagcaagaacaggagttcctggagaagaaatcagcaggagctgacggccgccaagaagcagcgcacccaaacaaagcgcggccctccaggcactgaagcgtagaagaggtatgagaagcagctggcgagatcgacggcacattatcaaccatcgagttccagcgggaggccctggagaatgccaaccaacacccgaggtgctcaagaacatgggtatgccgccaaggccatgaagcggcccatgacaacatggacatcgataaagttagtgagttatgcaggacattgctgaccagcaagaactgacagagagatttcaacagcaatttcaaacctgtagggtttggagaagagtttgacgaggatgagctcatggcgggaattagaagaactagaacaggaggaactagacaagaatttgcggaatcagtggaaccgaaacagtcctctaccaaatgttcctctatagccctaccatcaaaaccgccaagaagaagaagaggaggacgacgacatgaaggaaattggagaactggcctgatccatgtaa