

TMEM55A-transmembrane protein 55A Gene

Reference: PTXBC033892

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	type 2 phosphatidylinositol 4,5-bisphosphate 4-phosphatase; PtdIns-4,5-P(2) 4-phosphatase type II; ptdIns-4,5-P2 4-Ptase II; type 2 PtdIns-4,5-P2 4-Ptase; type II phosphatidylinositol 4,5-bisphosphate 4-phosphatase; transmembrane protein 55A		

For research use only.

GENE DESCRIPTION	transmembrane protein 55A
GENE ID	55529
ORIGIN SPECIES	Human
ACCESSIONS	BC033892
VECTOR	pDONR223
NCBI SYMBOL	TMEM55A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgctgatggggtggacgaacgctcgtctgtgtcagatcccaactcgggaatgtcactccaccgccccacgtacttgcagaagaagcagccccagagcggagctcccactccatatacagccattgccagtcagacgccagtggtattccagtaataaactgccgtgtgt
gccaatcactaatcaatttggatggcaagcttcaccagcatgtgttaagtgcacagtttgcaatgaagctacgccaatcaaaaccccccaacaggcagaataatgttagatgcccttgaattgtcttctcatttgaaggacacatctggcgaataggatgcccaagacccaactgtaga
cggataattaaccttggcccagtaatgttattctgaagaacaaccagctcagcctgcattgccaatccaaccagaaggtacaaaggtcgtgtgtgtggcactgtggaacacattcctgtggaatggaactgaggttcaacactctggcaaaatgccacactgcaaaaaatctcctcagtg
ggtagtgcacttcacgaagacgctcgtgtgtcatatattaccattggaatgatgtatttcattggagttgggttaactgttggcaccacagatttgcaggcgatttcgagcaacctatgtttcttggcgaattgcttctcctaggaattgatctgccttatccgagctgttattgggagccata
agagtcagttatccagaacacagttttgcatata