

### NOG-noggin Gene

Reference: PTXBC034027

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SYM1; SYNS1; SYNS1A; noggin; symphalangism 1 (proximal)		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	noggin
<b>GENE ID</b>	9241
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC034027
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	NOG
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

#### ORF (Only for reference, product provided without STOP codon)

```
atggagcgtccccagcctaggggtcacctctacgccctgggtggtcctggggctcggggcgaaccggccggcgccagcactatccacatccgcccggcaccagcagaacctgccccctggggacctatcgaacaccagacctatcttgacccaaggaaaag
gatctgaacgagacgtctgctgctcctcggggccactacgaccaggcttcatggccacctcgcctccgagaccggccggcggggggggggtgagctggggcgcgaggacctggcgagctggacagctctcggcagcggccgtcggggccatgc
cgagcgagatcaaaggctagatctccgagggctggcccaggcgaagaagcagcgcctaagcaagaagctcggaggaagtacagatgtgctgtgctgcagacattctccccgtgctgtacgctggaacgacctggcgagccgctttggccgctacgtgaaggtgg
gcagctctcagtaagcgtcgtcctccgtgcccaggcagtggtgtcaagccgtccaagctcctgacacctcaggtgctcgggtggcgctgcagcggcgcgggggccagcgtgctggctggattccatccaglaccccatcattccgagtgcaagtgctcgtcgtag
```