

KLHL3-kelch-like 3 (Drosophila) Gene

Reference: PTXBC034035

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PHA2D; kelch-like protein 3; kelch like family member 3		
For research use only.			

GENE DESCRIPTION	kelch-like 3 (Drosophila)
GENE ID	26249
ORIGIN SPECIES	Human
ACCESSIONS	BC034035
VECTOR	pDONR223
NCBI SYMBOL	KLHL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagggtgaaagtgtcaagctgagctccagactctgatacaggctgggatgatgagaagaaccagaggacgatcactgtcaaccctgccacatgggaaagcattcaaggttatgaatgaactcggagtaaacagctgtgtgtgacgtgatgtgtggcagaagatgtcgag
atagaagcccacccgtgtgtcctggcagcctgcagcccctactctgtgcgatgttcacaggtgacatgtctgagagtaaaagccaaaagatagaatcaaggacgtggatgggcagacgtctgagtaagctgattgactacatctactctgaaatcgaagtgactgaagagaatgtcca
gggtgctgtccggcagccagctgtgtcagctcatggatgttcggcagaactgtgtgacttctcagctcagtttgcatccacaaatgctgggcatccgtgcgtttgcagatgtacacacctgcactgacctctcagcagggccaatgcctacgcagagcagcacttccagaggt
gatgctaggagaagaatttcttagcctgagctgtgaccaggtgtgcagcttgatccagcgacaagctgaccgtttcttcagaagagaagggttttgaagctgtgatctcatggaatattgagaagaacccgttttagagcacatggcaagctgatggaacatgtccgacttctctctt
acctagggactacgtatgccaacggttgaagaagaagctttgataaagaataaacacacctgtaagacttctcattgagccatgaaataccatctctccctctggtacagagactattgattaagaaccaaggaccaagcccaggaactccagctcagcttccaagtag