



Reference: PTXBC034036

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BIRK-10; KCNJ13-PEN; KIR1.2; SESAME; ATP-sensitive inward rectifier potassium channel 10; ATP-dependent inwardly rectifying potassium channel Kir4.1; glial ATP-dependent inwardly rectifying potassium channel KIR4.1; inward rectifier K(+) channel Kir1.2; i		

For research use only.

GENE DESCRIPTION	potassium inwardly-rectifying channel, subfamily J, member 10
GENE ID	3766
ORIGIN SPECIES	Human
ACCESSIONS	BC034036
VECTOR	pDONR223
NCBI SYMBOL	KCNJ10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr