

LMO2-LIM domain only 2 (rhombotin-like 1) Gene

Reference: PTXBC034041

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RBTN2; RBTN1; RHOM2; TTG2; rhombotin-2; LIM domain only protein 2; LMO-2; T-cell translocation gene 2; T-cell translocation protein 2; cysteine-rich protein TTG-2; rhombotin-like 1; LIM domain only 2		

For research use only.

GENE DESCRIPTION	LIM domain only 2 (rhombotin-like 1)
GENE ID	4005
ORIGIN SPECIES	Human
ACCESSIONS	BC034041
VECTOR	pDONR223
NCBI SYMBOL	LMO2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctcggccatgaaaggaagagcctggaccctcagaggaaaccagtggatgaggtgctgcagatcccccatccctgctgacatgcggcgctgccagcagaacattggggaccgctacttctgaaggccatcgaccagtactggcacgaggactgctgagctgcgacctctgtggctgccggctgggtgaggtggggcggcgcctctactacaaactggccggaagctctgccggagagactatctcaggcttttgggcaagacgggtctctgcgcacccctgtgacaagcgggattctgtcctatgagatgacaatgcgggtgaaagacaaagtgtatcacctggaatgtttcaagtgcggcctgtcagaagcattctgttaggtgacagatacctctcatcaactctgacatagttgtgcgaacaggacatctacgagtggaactaatggaatgatag