



Reference: PTXBC034047

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FIGLER7; LRANK1; LRRN5; leucine-rich repeat neuronal protein 2; fibronectin type III, immunoglobulin and leucine rich repeat domain 7; glioma amplified on chromosome 1 protein; leucine rich and ankyrin repeats 1; leucine rich repeat neuronal 5; leucine-ri		

For research use only.

GENE DESCRIPTION	leucine rich repeat neuronal 2
GENE ID	10446
ORIGIN SPECIES	Human
ACCESSIONS	BC034047
VECTOR	pDONR223
NCBI SYMBOL	LRRN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]