

DNAJC27-DnaJ (Hsp40) homolog, subfamily C, member 27 Gene

Reference: PTXBC034049

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RBJ; RabJS; dnaJ homolog subfamily C member 27; DnaJ (Hsp40) homolog, subfamily C, member 27; Ras-associated protein Rap1; rab and DnaJ domain-containing protein; DnaJ heat shock protein family (Hsp40) member C27		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 27
GENE ID	51277
ORIGIN SPECIES	Human
ACCESSIONS	BC034049
VECTOR	pDONR223
NCBI SYMBOL	DNAJC27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggccaacatgccgaagcgaagcctggcaggctctccgcacaaagtcctccatgggcaacgccgaagtgggaaaaagctgtattataaagcgatactgtgagaaaagattcgtctctaaatcctggcaacaattggaattgactatggagtcacaaaggtacacgtca
gagacagagaaatcaaatgtaacatcttggatggcctggacatccctctctatgaggttcgaaatgagtttacaaggacacacaggggtgatactggctctatgatgtgggcagaaagactccttggaccccttgatgcgtggcctggcagaaatgaagcaagagctggacatcagga
aacatggaaaatattatattttagtttgccaacaaagattgattgtacaaacatcgtgtgtagatgaaagtgaaggacgtcttggcctgaaagcaaaaggggttcctgtacttgaacttcagcacaaactggagaaggcattaatgagatgttcagacacctttatataatccatagtgatttat
gtgaaaatggcgggaaacgccctaccaccaatagcagtgctagtttaccacaagaacagcagatgccattcgagaattcgaaatagtaaacagagttgggacatcgtgggagtcacaaactggggcctcaagggatgaagtcaataaagcgtatcggaaacttgctgtcttccccc
tgacaaatgtgtagcctggcagtgagatgccctcaaaagcagttgtgaatgctcgacagccctcctgaaaaacatcaagtag