

TM4SF1-transmembrane 4 L six family member 1 Gene

Reference: PTXBC034145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H-L6; M3S1; TAAL6; transmembrane 4 L6 family member 1; membrane component, chromosome 3, surface marker 1; transmembrane 4 superfamily member 1; tumor-associated antigen L6; transmembrane 4 L six family member 1		

For research use only.

GENE DESCRIPTION	transmembrane 4 L six family member 1
GENE ID	4071
ORIGIN SPECIES	Human
ACCESSIONS	BC034145
VECTOR	pDONR223
NCBI SYMBOL	TM4SF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgctatgggaagtgtgcacgatgcacgtgacattctctgttgggtcgcgcctcctgtgcatcgcggctaataatttgccttacttcccaatggggaacaaagtatgcctccgaaaccacctcagccgttcgtgtgttcttctggcatcgtaggaggtggcctgctgatgctcctgc
cagcatttgcctcattgggtggaacaggatgactgctgtgggtgctgtggccatgaaactgtggcaaacgatgtgcgatcttctctgtattggctgctctcattggaattgcaggatcgtggctactgtgcatgtgtggcagccctggccttagcagaaggaccactatgcttgaattccctcg
gccagtggaaactacaccttggcagcaccgagggccagtacctctggatacctccacatgtgccgagtgcatgaaccaagcacattgtggaatggaatgtatctctgtttctatcctcttggctcttgggtgaattgaattcatctgtgtcttattcaagtaataaatggagtgcttggaggca
tatgtggcttttgcctctaccaacagcaaatgactgctaaagaaccaaccagagacagagccacaatcttctctatttcattgtaatttatatttcttgcatttcattgtataaacttgcattag