

MGC39821-hypothetical protein MGC39821 Gene

Reference: PTXBC034236

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	long intergenic non-protein coding RNA 663		

For research use only.

GENE DESCRIPTION	hypothetical protein MGC39821
GENE ID	284440
ORIGIN SPECIES	Human
ACCESSIONS	BC034236
VECTOR	pDONR223
NCBI SYMBOL	MGC39821
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtcaccagaggctcaagcccttctcaccatggtaccacctctcgaagctgcctccaaagccctctagctgtgtgcagccacaacaagtgcaccaccaggggcttccagctgcaggtgcaggttccagagccgggactggggaacctgcagtcacacctcagc
aattctggtggcagccagccagccagggcagcgtgatgaccgtaataaaaagtgggatgacagctccacactgggtgatgagttgagaacaccagggcagaggtacgtggtcttggatggttcagtgagactgcactacaattgcagtgcaaaagtgcagagttga