

TSPAN9-tetraspanin 9 Gene

Reference: PTXBC034280

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET-5; NET5; PP1057; tetraspanin-9; new EST tetraspan 5; tetraspan NET-5; transmembrane 4 superfamily member tetraspan NET-5; tspan-9; tetraspanin 9		

For research use only.

GENE DESCRIPTION	tetraspanin 9
GENE ID	10867
ORIGIN SPECIES	Human
ACCESSIONS	BC034280
VECTOR	pDONR223
NCBI SYMBOL	TSPAN9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggggaggccacatttgagcaagctgccctgtccctggcctgtgtaaggagtttcagccgctccgcaggaaccccaagtcagattccggagcagacacatccggcggagagactcagcagacaagtgcagttgcacgggtggaccggg
gcctcgtgcgtttttctgtgtgggtgggtgggtgtgtttatgcctatcaatgcaatttttaattttgtaataatcaacagcaaaagcctagtgcatfgggagatgtgcaacctccctgaaaatctttctgtttctggagtactcagggtggcctctgccccagagcctttgccagtgct
cccaccagccccacatcctcgtctgtttgcagagcctcatctacaggtcccacgtgcctctttactcactctgcgcttggcgtttgtatttggcttag