

UBE2Q2-ubiquitin-conjugating enzyme E2Q family member 2 Gene

Reference: PTXBC034342

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ubiquitin-conjugating enzyme E2 Q2; E2 ubiquitin-conjugating enzyme Q2; ubiquitin carrier protein Q2; ubiquitin conjugating enzyme E2Q family member 2; ubiquitin-protein ligase Q2; ubiquitin conjugating enzyme E2 Q2		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2Q family member 2
GENE ID	92912
ORIGIN SPECIES	Human
ACCESSIONS	BC034342
VECTOR	pDONR223
NCBI SYMBOL	UBE2Q2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaactctcacaacagggctggagcagtgccctactcaatagaatcggatcatgcaataaatgccaccttagtcaaaggcaaagccagagtgagttggagcaataagaatcaatataatctagcaagagcccaacaatcctataatccattgtacagatacatgagaaaaat
ggctgtacacccctcaaaaggaatggctaa