

YSK4-YSK4 Sps1/Ste20-related kinase homolog (S. cerevisiae) Gene

Reference: PTXBC034417

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YSK4 Sps1/Ste20-related kinase homolog; SPS1/STE20-related protein kinase YSK4; YSK4; RCK; mitogen-activated protein kinase kinase kinase 19; regulated in COPD kinase; regulated in COPD, protein kinase; yeast Sps1/Ste20-related kinase 4		

For research use only.

GENE DESCRIPTION	YSK4 Sps1/Ste20-related kinase homolog (S. cerevisiae)
GENE ID	80122
ORIGIN SPECIES	Human
ACCESSIONS	BC034417
VECTOR	pDONR223
NCBI SYMBOL	YSK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttgttctggtggctcaatcttagtattataaacgtttggccattgctgagatggtgttctgtaaatatacgaacaaatactcaagggtgttcttatccatgagaactgtgtgtacatcgcatatcaaaaggaaataatgttatgctcatgccaactggaataaaagctgat
tgactttggctgtgccaggcgtttggcctgggcaggtttaaatggcaccacagtgacatgcttaagtcctatgcatgggaactccatattggatggcccagaatgcatcaatgagctggctatggacggaaatcagatactggagcattggtgtactgtgtttgagatggctacaggggaagc
ctccactggcttccatggacaggatggccgcatgtttatcggagcacaccgagggctgatgctcctttaccagaccacttctcagaaaaatgcagcagactttgtgcgcatgtgcctgaccaggttaa