



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FTMT-ferritin mitochondrial Gene

Reference: PTXBC034419

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MTF; ferritin, mitochondrial; ferritin H subunit; mitochondrial ferritin; ferritin mitochondrial		

For research use only.

GENE DESCRIPTION	ferritin mitochondrial
GENE ID	94033
ORIGIN SPECIES	Human
ACCESSIONS	BC034419
VECTOR	pDONR223
NCBI SYMBOL	FTMT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1