

C1orf190-chromosome 1 open reading frame 190 Gene

Reference: PTXBC034422

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NF-kappa-B activator C1orf190; C1orf190; LRAP35a; LRP35A; leucine rich adaptor protein 1; leucine repeat adapter protein 35A; leucine repeat adaptor protein 35a		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 190
GENE ID	541468
ORIGIN SPECIES	Human
ACCESSIONS	BC034422
VECTOR	pDONR223
NCBI SYMBOL	C1orf190
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggggaccgtggagctccagacgcctgacctgcgggatgtggagggttaagtgaggcagggaagaccctgaagggtctccgagggtgcgaggcagtggtgagctgggaacctctggcgcctctgctccaggggcgtctagcaccggccacgacttggggacaag
atcatggcgtgaagatggagctggcttacctgcgagccatcgaatgaagatcctgcagcagctggtgacctgaatgagggtcatcgaggcagtgctgctgttggaggagcggggagcttgaccagtcattgcagcagcctcaccagcagtcattatagcctgacaggcggga
gcccaggccgtcaaggcagggcagctgggacagcctgccagaccagcaccgaccggtggacagtgctctatggcagcttcctggacacagtggtggccagcagctggatgaacaggggccacctggggctccagctccagatggactgggcaaaagttagct
ggtggagagagggccaggactgaggtgatgtggcagccaccaggtaggagcttgagagctgtgtggaagccccaggggagaggtcgaaggtggaccacctgagtcaccagaggatgagagtccaagctgggttcgagggccactgggtctggagcagtgccaggat
gatgtgaccttctgtaa