



Reference: PTXBC034479

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERCC6-PGBD3 fusion protein; ARMD5; CKN2; COFS; COFS1; CSB; POF11; RAD26; UVSS1; DNA excision repair protein ERCC-6; ATP-dependent helicase ERCC6; Cockayne syndrome group B protein; cockayne syndrome protein CSB; excision repair cross-complementation group		

For research use only.

GENE DESCRIPTION	piggyBac transposable element derived 3
GENE ID	267004
ORIGIN SPECIES	Human
ACCESSIONS	BC034479
VECTOR	pDONR223
NCBI SYMBOL	PGBD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]