

ZNF738-zinc finger protein 738 Gene

Reference: PTXBC034499

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			

For research use only.

GENE DESCRIPTION	zinc finger protein 738
GENE ID	148203
ORIGIN SPECIES	Human
ACCESSIONS	BC034499
VECTOR	pDONR223
NCBI SYMBOL	ZNF738
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagacctgaggtatggagtgatcctgtcaagggggcaagtggataccctgggctgagaggaatctctggagtactctatftgaaaagggcctgacatttagggatgtgcatagaattctctcaggaggagtgcaatgcctggacactgctcagcaagattgtatagga
aagtgatgttagagaactcagaaacctgggtctctgggtattgatgtcttaagccagatctgatcactctctggagcaagggaaaagatcctggaaatagaagagacacagtaggtagccacacccccagaaaggggtattgaaatttctactataattatattgctgtcaggtgttc
ttccagctctgcaatattgctttatataattggaaccttag