



Reference: PTXBC034502

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	choline dehydrogenase, mitochondrial; CDH; CHD; choline dehydrogenase		
For research use only.			

GENE DESCRIPTION	choline dehydrogenase
GENE ID	55349
ORIGIN SPECIES	Human
ACCESSIONS	BC034502
VECTOR	pDONR223
NCBI SYMBOL	CHDH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]