

KRTAP3-2-keratin associated protein 3-2 Gene

Reference: PTXBC034582

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP3.2; KRTAP3.2; keratin-associated protein 3-2; high sulfur keratin-associated protein 3.2; keratin associated protein 3.2; keratin associated protein 3-2		

For research use only.

GENE DESCRIPTION	keratin associated protein 3-2
GENE ID	83897
ORIGIN SPECIES	Human
ACCESSIONS	BC034582
VECTOR	pDONR223
NCBI SYMBOL	KRTAP3-2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattgctgtgcctctcgcagctgcagtgctcccaactgggctgccaccacatctgtctcgcacaaatctgccgtgtggagtgctgtgccagcactgccacacacagtttggtactggagccatctgctgtgacaactgtccccaccctgccacattctcagccctgcg
tgccccactgtctctcgaactctgcagccaactccggcctggagaccctcaacctcaccacttcaactcagccctgctgtgagccctgctcccaagaggtgctga