



Reference: PTXBC034607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ASP; Calmbp1; MCPH5; abnormal spindle-like microcephaly-associated protein; asp (abnormal spindle) homolog, microcephaly associated; abnormal spindle microtubule assembly		

For research use only.

GENE DESCRIPTION	asp (abnormal spindle) homolog, microcephaly associated (Drosophila)
GENE ID	259266
ORIGIN SPECIES	Human
ACCESSIONS	BC034607
VECTOR	pDONR223
NCBI SYMBOL	ASPM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]