

## LYPD4-LY6/PLAUR domain containing 4 Gene

Reference: PTXBC034629

|                |                                                                                                                                      |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | SMR; ly6/PLAUR domain-containing protein 4; sperm membrane receptor; testicular tissue protein Li 159; LY6/PLAUR domain containing 4 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | LY6/PLAUR domain containing 4                                                                     |
| <b>GENE ID</b>          | 147719                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC034629                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | LYPD4                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgggacccagcattgagacttgcagctgttgccttctagggccatctccactctgctctgtatgtctgtgggctggatgctataagaccagaaaggactgcaagggagtcgtgggtttaaggctgcagctcgtctcttaccctgcgcaaatctctacattgttcc  
ccacccggagtgctcctacagtcgctctgctgcttctctgcaacaacctcaccaattggagcctttgtgaaactcaaggccagcactcctaagctatcacatctgcctctgtagctgcccgaacctgtgtggcgagcacatgaaggattgctcccaattttgcacc  
actaatcttgccttggctgctctacgtgttacagttccaccttaaaattcaggcaggtttctcaataccacctctcctcatggggtgtgctcgtgaacataaccagcttttagcagattttcatcatattgggagcatcaagtgactgaggtcctcaacatcttagagaagtcagattgtt  
ggtgcagcatctccaggcaagatctgcttgggtgtcgtcttaggcctctgtttgcctcagggaactga