

CDK2AP1-cyclin-dependent kinase 2 associated protein 1 Gene

Reference: PTXBC034717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DOC1; DORC1; ST19; doc-1; p12DOC-1; cyclin-dependent kinase 2-associated protein 1; CDK2-associated protein 1; Deleted in oral cancer-1; deleted in oral cancer 1; cyclin dependent kinase 2 associated protein 1		

For research use only.

GENE DESCRIPTION	cyclin-dependent kinase 2 associated protein 1
GENE ID	8099
ORIGIN SPECIES	Human
ACCESSIONS	BC034717
VECTOR	pDONR223
NCBI SYMBOL	CDK2AP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcttacaaccgaactggccgcacatgccgcgcgcctcaacgcgcgtgggagtggtccactgccttcaccagcatggcaacgtcttcacagtaccgccagctgctcagtgactacggccaccgtccctaggtacaccaggggaactgggaacagccaggtgcccc
aaagcaaatacgcggagctgctggccatcattgaagagctggggaaggagatcagaccacgtacgcaggagcaagagtgccatggagagctgaagcgcggcatcattcacgtagaggactggttcgggagtgcttggcagaacggaacggaatgccagatcctag