

## EPB41L2-erythrocyte membrane protein band 4.1-like 2 Gene

Reference: PTXBC034718

|                |                                                                                                                                                                    |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | 4.1G; band 4.1-like protein 2; erythrocyte membrane protein band 4.1 like-protein 2; generally expressed protein 4.1; erythrocyte membrane protein band 4.1 like 2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | erythrocyte membrane protein band 4.1-like 2                                                     |
| <b>GENE ID</b>          | 2037                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC034718                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | EPB41L2                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgactactgaagtaggctctgtgtgaaggaaggactctagccagttggaacagatgcaaccaaggaaaaactaaagaagtagcagaaatcagcagaatcagctctccgatccagaggaggaaaaaggttccagccacctctcagctgaaagccaaagtagtctacg  
ccgccagaagagagagaaggaaacatcgagagcagggtgtttctcgggtcattaccgcatggcttaagaagcaaaagtcatacttagtagtgccaaagatggaggagataaaaaagcctacccaagctgtgttgagaacagctcttagataaagaggaaacccctccaga  
agaacagagacaggctaaagggtgatgctgaagaatggctcagaagaacaagagattaaagttgaagtcagggaagaaaaacccctcagtgagcaagggaagaaacccctcagtgagcaaggtggagatgcagccctactgaattagtaagtaggagagagaaggttaaag  
aaacacagggaagacaattagaaggaggagcagcaaaaagggaaccaaggagtcagacacaaatgagctgaagcagagaaagcagctcaaaaagtcaccaagaagaccaaaactgtccagtgtaaaagtgaacctcttagatggcaccgaatacagctgtgacctggagaacat  
gccaaggagacaagtgatttgaacaagtgtgtgaacacctcaatctctggagaaagactacttggactttgtttcaggaaagccctgagcagaaaaacttggttagatcctgttaagaaataaagagacaactgagaacacctccatggctatctactttaaagtttatcctctgat  
ccttctcaattgactgaagatacaccagatacttctgtgcttcagctccggcaggacattgcctctggccgctccctgtcttctgtgactcatgctctcctgggatcctacacccctgcaggctgaacttggtgactatgaccagaagaacatggcagcatcgacctcagtgaaatccagt  
ttgcccctactcagactaaggagctggagagaaggtggcagagctgcacaaaacccacaggggcttatcgccagcacaagctgattccagcttcttagaaaaatgcaaaagaggctttccatgtatggtgtgtgacatcatcatccaaggaactcagaaggtgtgtgacatcaagctggcggt  
gtgtgtaattgacttctcattacaagacagactgcgaalcaatcgtttgtgtggccgaaactttaaatttctataaacgcagtaactctacataaagtcagacggcagagctggaaagtttgagagtaccattgattcaaaactgccaaccaccggggcagcgaagaaactatg  
gaaagtgtgctggagcatcatacttctacaggctgtttctcagagcagccacaaaagccaagttcctgacctgggggtccaaatttcgctatagtgccgcacccaagcacagaccgccaaggccagcaccctcatagataggccagcaccactttgagcgactctagtaaac  
gggtctccaggagctctagatggagctccgattgggtgtcatggaccaaaagttctatgaaggatttctcgtgcgtgtgggagatttcagcctatggacctggactgtcagcattgcccgtgtgtacaagatggggacggcagaggaggagtagagaagcccaactaaagccccacatttga  
gctcattgaaggaaagagttcagactctgaatatagtggagagaagagcgggcagaggttgggaaagacgaaagagtaa