



Reference: PTXBC034737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MPDZ; PDZRN2; E3 ubiquitin-protein ligase LNX; PDZ domain-containing ring finger protein 2; ligand of numb-protein X 1, E3 ubiquitin protein ligase; multi-PDZ-domain-containing protein, E3 ubiquitin-protein ligase LNX; numb-binding protein 1; ligand of nu		

GENE DESCRIPTION	ligand of numb-protein X 1
GENE ID	84708
ORIGIN SPECIES	Human
ACCESSIONS	BC034737
VECTOR	pDONR223
NCBI SYMBOL	LNK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]