



Reference: PTXBC034774

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C6orf10; TSBP; testis specific basic protein; chromosome 6 open reading frame 10		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 10
GENE ID	10665
ORIGIN SPECIES	Human
ACCESSIONS	BC034774
VECTOR	pDONR223
NCBI SYMBOL	C6orf10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]