

ZDHC19-zinc finger, DHHC-type containing 19 Gene

Reference: PTXBC034801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DHHC19; DHHC-19; zinc finger DHHC domain-containing protein 19; zinc finger, DHHC domain containing 19; zinc finger DHHC-type containing 19		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 19
GENE ID	131540
ORIGIN SPECIES	Human
ACCESSIONS	BC034801
VECTOR	pDONR223
NCBI SYMBOL	ZDHC19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacactttaacgatgccagccgctggtgaaggagccccatccccctgctctggtccacgtccctggtctccctagcctcttctgctccctcaatggtgctgctggtctcttttcagtggtcctcttcttcgcatccctctgtttcactcaactctcagaccctggcatcttacatcaagg
ctccgctgagcaggcccccttgacggtgacagtggtgtggtgaaccacgggctccgctgcaatggtgtccaaagtgtgcttccaccgccccggactaccactgccccctggtgcaacatctgtgtggaggtgagtgccccctacctgtccaccagccatgccccacat
ctggctgggaccaggagcctccctacacaggggcagggtgggtgtgaagctgtttcaggactgctaaaggccaagggtggcagtaactccagccacccttggtcactcatgtgttctccactgcacctgtctctgggaggctggccgttcctgcagaggtgctgacatgcacag
ggcatccgcatctgtctggtggtcaggccaggggaggagggcagccctctgtccacccctgggtgctgaacacgagataaccaagccccctcaggctggctgcccgcctggtctgggaggccttgctgtctcgcctgggtgtccaagttggagggtctctggtctgagt
ga