



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

EDA2R-ectodysplasin A2 receptor Gene

Reference: PTXBC034919

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

EDA-A2R; EDAA2R; TNFRSF27; XEDAR; tumor necrosis factor receptor

superfamily member 27; EDA-A2 receptor; X-linked ectodysplasin-A2 receptor;

ectodysplasin A2 isoform receptor; tumor necrosis factor receptor superfamily member

XEDAR; ectodysplasin A2 recep

For research use only.

GENE DESCRIPTION	ectodysplasin A2 receptor
GENE ID	60401
ORIGIN SPECIES	Human
ACCESSIONS	BC034919
VECTOR	pDONR223
NCBI SYMBOL	EDA2R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagattgccaaaataatgagtactgggccaatggggcgggtgtgtctacctgccaacgggtgtgctctggacaggagctatcaagaattgtgtgtatggagagggtggagatgctactgcacagcctgcctctctgcaggctacaaaagcagctggggccacacagatgtcagagtgcacatcctgtgctgtcatcaatctgtgtcagaaggtcaactgcacagctacctctaattgctgtctgtggggactgtttgccagggttaccgaagaacacgcattggaggcctcaggagcaagaagtgcattccctgcacgaagcagacccccacctctgaggttcaatgtgccttcagttgagcttagtggagcagatgagctggcctccccctcaggagcccaactgttgacgttggagcagctctgactgtgtttactctgcgcctctcctgcctctcttctctactcaagaacgattcttcaacagagctcagcagctggtagctgttggctgagcttggagctgtataaacagcaaaaagggaatctcttccccctgcaccaggcaggagcagctgtgagtgtgccaaagtgatgagcaactcttcagacagccactctcggaggcagctgcagctgactgtgactgtgtctccccacagaggtcttcttactgtcctctgcacctcagagagccactccctctgcacagccccctgaatgcacagagctggcagctccaaaattttccagctctctctctatactgagctgcagactctgtggggaaacacagctgaagacgctgaagacagctggagcagctggagctcaatgtcttgaatttccacgctctaa

More info: Please contact us for any question: contact@proteogenix.fr